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


- [1] Andrew C Faust , Pearl Rajan, Lyndsay A Sheperd, Carlos A Alvarez, Phyllis McCorstin, Rebecca L Doebele Impact of an Analgesia-Based Sedation Protocol on Mechanically Ventilated Patients in a Medical Intensive Care Unit *Anesth Analg*, 903-9 Oct 2016
- [2] B V Christensen , L P Thunedborg Use of Sedatives, Analgesics and Neuromuscular Blocking Agents in Danish ICUs 1996/97. A National Survey *Intensive Care Med* 25 (2), 186-91 Feb 1999
- [3] European Heart Journal (1996) 17, 354-381 Guidelines Heart rate variability Standards of measurement, physiological interpretation, and clinical use Task Force of The European Society of Cardiology and The North American Society of Pacing and Electrophysiology (Membership of the Task Force listed in the Appendix)
- [4] Jeanne M, Logier R, De Jonckheere J, Tavernier B. Validation of a graphic measurement of heart rate variability to assess analgesia/nociception balance during general anesthesia, *IEEE Proceedings* 2009.
- [5] Jeanne M, Delecroix M, De Jonckheere J, Keribedj A, Logier R, Tavernier B. Variations of the Analgesia Nociception Index during propofol anesthesia for total knee Replacement. *Clin J Pain* 2014 Dec; 30(12):1084-8.
- [6] Funcke et al. Validation of Innovative Techniques for Monitoring Nociception during General Anesthesia. *Anesthesiology* 2017 Aug;127(2):272-283
- [7] Logier R, et al. Heart rate variability analysis for arterial hypertension etiological diagnosis during surgical procedures under tourniquet. *Conf Proc IEEE Eng Med Biol Soc* 2011;2011:3776-9.
- [8] Daccache G, et al. A targeted remifentanil administration protocol based on the analgesia nociception index during vascular surgery. *Anaesth Crit Care Pain Med* 2017 Aug;36(4):229-232.
- [9] Boselli E, et al. Prediction of immediate postoperative pain using the analgesia/nociception index: a prospective observational study. *Br J Anaesth* 2014 Apr;112(4):715-21.
- [10] Broucqsault-Dédrie C, De Jonckheere J, Jeanne M, Nseir S. Measurement of Heart Rate Variability to assess pain in sedated critically ill patients: A prospective observational study. *Plos One* 2016; 11(1): e0147720
- [11] Boselli E, Fatah A, Ledochowski S, Allaouchiche B. ANI and BIS variations in supine and prone position during closed-tracheal suction in sedated and myorelaxed ICU patients with severe COVID-19: A retrospective study. *J Clin Monit Comput*. 2020 Nov 6 : 1-7.
- [12] Ali Jendoubi, Ahmed Abbes, Salma Ghedira, and Mohamed Houissa. Pain Measurement in Mechanically Ventilated Patients with Traumatic Brain Injury: Behavioral Pain Tools Versus Analgesia Nociception Index *Indian J Crit Care Med*. 2017 Sep; 21(9): 585-588.
- [13] Abdullayev R, Uludag O, Celik B. Analgesia Nociception Index: assessment of acute postoperative pain. *Rev Bras Anesthesiol*. 2019 Jul - Aug;69(4):396-402. doi: 10.1016/j.bjan.2019.01.003. Epub 2019 Aug 6.
- [14] Devlin JW et al., Executive Summary: Clinical Practice Guidelines for the Prevention and Management of Pain, Agitation/Sedation, Delirium, Immobility, and Sleep Disruption in Adult Patients in the ICU. *Crit Care Med*. 2018 Sep;46(9):1532-1548.
- [15] Carlos Chamorro, Miguel Angel Romera, Sara Alcantara. ICU Department Hospital Universitario Puerta de Hierro. Protocolo de sedoanalgesia para los pacientes con COVID19.
- [16] <https://semicyuc.org/wp-content/uploads/2020/03/PROTOCOLO-ANALGOSED-COVID-SEMICYUC.pdf>
- [17] Henry D. Upton, MBBS, BMedSc (Hons), Guy L. Ludbrook, MBBS, FANZCA, PhD, Andrew Wing, BMBS (Hons), BSc (Hons), FANZCA, and Jamie W. Sleight, MD Intraoperative "Analgesia Nociception Index"- Guided Fentanyl Administration During Sevoflurane Anesthesia in Lumbar Discectomy and Laminectomy: A Randomized Clinical Trial *Anesthesia & Analgesia* July 2017 doi: 10.1213/ANE.0000000000001984

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ANI Monitor V1, ANI Monitor V2 y ANI MOC-9 son dispositivos médicos de clase IIa, fabricados por Mdoloris Medical Systems y han sido evaluados para la CE por Bureau Veritas Italia (1370). ANI Monitor V1, ANI Monitor V2 y el nombre y el logotipo de ANI MOC-9, ANI MR son marcas comerciales registradas. © 2021 MDoloris Medical Systems. Reservados todos los derechos. MD/QUA/EN14.2 (es) v.01
MDoloris Medical Systems SAS